

इंटरनेट

मानक

Disclosure to Promote the Right To Information

Whereas the Parliament of India has set out to provide a practical regime of right to information for citizens to secure access to information under the control of public authorities, in order to promote transparency and accountability in the working of every public authority, and whereas the attached publication of the Bureau of Indian Standards is of particular interest to the public, particularly disadvantaged communities and those engaged in the pursuit of education and knowledge, the attached public safety standard is made available to promote the timely dissemination of this information in an accurate manner to the public.

“जानने का अधिकार, जीने का अधिकार”

Mazdoor Kisan Shakti Sangathan

“The Right to Information, The Right to Live”

“पुराने को छोड़ नये के तरफ”

Jawaharlal Nehru

“Step Out From the Old to the New”

IS 5952 (1970): Methods for determination of parathion residues [FAD 1: Pesticides and Pesticides Residue Analysis]



“ज्ञान से एक नये भारत का निर्माण”

Satyanarayan Gangaram Pitroda

“Invent a New India Using Knowledge”



“ज्ञान एक ऐसा खजाना है जो कभी चुराया नहीं जा सकता है”

Bhartrhari—Nitiśatakam

“Knowledge is such a treasure which cannot be stolen”

BLANK PAGE




IS : 5952 - 1970

Indian Standard

METHOD FOR DETERMINATION
OF PARATHION RESIDUES

UDC 633.1:632.95.028:543.432



© Copyright 1971

INDIAN STANDARDS INSTITUTION
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG
NEW DELHI 1

November 1971

Gr 2

*Indian Standard***METHOD FOR DETERMINATION
OF PARATHION RESIDUES**

Pest Control Sectional Committee, AFCDC 6

Chairman

DR K. D. PAHARIA

Representing

Directorate of Plant Protection, Quarantine & Storage (Ministry of Food, Agriculture, Community Development & Co-operation)

Members

SHRI M. L. ADVANI

DR K. KALYAN (*Alternate*)

DR P. N. AGARWAL

DR K. BAGCHI

DR H. L. BAMI

DR N. K. BASU

DR T. GUHA (*Alternate*)

MAJ-GEN M. S. BOPARAI

DIRECTOR OF AGRICULTURE

DR P. L. CHATURVEDI (*Alternate*)

DIRECTOR OF STORAGE & INSPECTION

DEPUTY DIRECTOR (STORAGE, TECHNICAL) (*Alternate*)

DRUGS CONTROLLER

SHRI G. D. GOKHALE

SHRI V. V. KETKAR (*Alternate*)

DR R. GONCALVES

SHRI H. J. DADACHANJI (*Alternate*)

DR K. C. GULATI

SHRI T. GEORGE JOHN

DR S. N. NAG (*Alternate*)

DR LALLAN RAI

SHRI M. K. ZUTSHI (*Alternate*)

SHRI S. K. MAJUMDAR

SHRI M. MUTHU (*Alternate*)

Tata Fison Industries Ltd, Bombay

Ministry of Defence (R & D)

Ministry of Health, Family Planning, Works, Housing & Urban Development

Central Forensic Science Laboratory, Central Bureau of Investigations, New Delhi

I.C.I. (India) Private Ltd, Calcutta

Directorate General of Armed Forces Medical Services, Ministry of Defence

Department of Agriculture, Government of Uttar Pradesh

Directorate of Storage & Inspection (Ministry of Food, Agriculture, Community Development & Co-operation)

Indian Pharmacopoeia Committee (Ministry of Health, Family Planning, Works, Housing & Urban Development)

Bombay Chemicals Private Ltd, Bombay

National Organic Chemical Industries Limited, Bombay

Indian Agricultural Research Institute, New Delhi

Hindustan Insecticides Limited, New Delhi

Cyanamid India Limited, Bombay

Central Food Technological Research Institute (CSIR), Mysore

(*Continued on page 2*)

INDIAN STANDARDS INSTITUTION
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG
NEW DELHI 1

(Continued from page 1)

<i>Members</i>	<i>Representing</i>
SHRI L. P. MATHUR	Ministry of Defence (DGI)
SHRI K. P. SINGH (<i>Alternate</i>)	
DR P. R. MEHTA	Pesticides Association of India, New Delhi
SHRI K. S. MEHTA	Bharat Pulverising Mills Pvt Ltd, Bombay
SHRI S. CHATTERJI (<i>Alternate</i>)	
SHRI D. N. NAKHATE	Bayer (India) Ltd, Bombay
DR D. V. KASHELIKAR (<i>Alternate</i>)	
SHRI R. NARASIMHAN	Directorate General of Technical Development (Ministry of Industrial Development and Internal Trade)
SHRI S. L. PATEL	American Spring & Pressing Works Private Ltd, Bombay
SHRI L. M. PATEL (<i>Alternate</i>)	
PLANT PROTECTION OFFICER	Department of Agriculture, Government of Mysore
DR S. PRADHAN	Indian Agricultural Research Institute, New Delhi
DR RATTAN LAL (<i>Alternate</i>)	
SHRI P. RAMDAS	Central Warehousing Corporation, New Delhi
DR A. SARUP	Pest Control Corporation, Calcutta
DR VINOD R. SHAH	Pest Control (India) Pvt Ltd, Bombay
SHRI L. R. SHARMA	Shaw Wallace & Company Limited, Calcutta
SHRI N. R. JAYASUNDAR (<i>Alternate</i>)	
DR K. S. SINGH	Indian Veterinary Research Institute, Izatnagar
DR R. P. CHAUDHURI (<i>Alternate</i>)	
DR B. L. WATTAL	National Institute of Communicable Diseases, Delhi
DR HARI BHAGWAN, Deputy Director (Agri & Food)	Director General, ISI (<i>Ex-officio Member</i>)

Secretary

SHRI E. N. SUNDAR
Deputy Director (Agri & Food), ISI

Pesticide Residue Analysis Subcommittee, AFCDC 6 : 7

Convener

DR H. L. BAMI
Central Forensic Science Laboratory, Central Bureau
of Investigations, New Delhi

Members

DR K. BAGCHI	Central Committee for Food Standards, New Delhi
SHRI D. S. CHADHA (<i>Alternate</i>)	
SHRI K. BALASUBRAHMANIAN	King Institute, Madras
DR S. L. CHOPRA	Punjab Agricultural University, Ludhiana
SHRI H. J. DADACHANJI	National Organic Chemical Industries Limited, Bombay
SHRI P. CHAWLA (<i>Alternate</i>)	
DR R. S. DEWAN	Indian Agricultural Research Institute, New Delhi
SHRI V. V. KETKAR	Bombay Chemicals Private Ltd, Bombay
DR K. KRISHNAMURTHY	Directorate of Storage & Inspection (Ministry of Food, Agriculture, Community Development & Co-operation)

(Continued on page 8)

Indian Standard

METHOD FOR DETERMINATION OF PARATHION RESIDUES

0. FOREWORD

0.1 This Indian Standard was adopted by the Indian Standards Institution on 29 December 1970, after the draft finalized by the Pest Control Sectional Committee had been approved by the Agricultural and Food Products Division Council and the Chemical Division Council.

0.2 Due to the increased use of pesticides in the country, the question of their harmful residues and their effect on public health has been engaging the attention of the people concerned with public health. With a view to helping in collection of authentic data on pesticidal residues and in order to enable the health authorities to implement the existing Prevention of Food Adulteration Rules relating to the limits of pesticidal residues, the Indian Standards Institution has undertaken the task of preparing a series of Indian Standard methods of test for estimation of pesticidal residues.

0.3 While preparing these standards, due attention has been paid to the methods of test being followed in various laboratories of the country. For the sake of wider adoption, only those methods have been prescribed in these Indian Standards which, as far as possible, avoid the use of costly equipment and reagents.

0.4 Prior to the estimation of pesticidal residues, it is necessary to identify qualitatively the pesticidal residues. Also, it is important that the solvents that are used in these methods of test are pure and solvents are purified before analysis. Taking into consideration these two important factors, two standards covering these fields are under preparation. These standards would be necessary adjuncts to this standard.

0.5 Parathion formulations have become popular in the agricultural and public health field for the control of various pests and the method of test for estimating its residues as described in the standard is one in the series of the Indian Standard methods of test for estimating residues of various pesticides in different foodstuffs.

0.6 In reporting the result of a test or analysis made in accordance with this standard, if the final value, observed or calculated, is to be rounded off, it shall be done in accordance with IS:2-1960*.

*Rules for rounding off numerical values (*revised*).

1. SCOPE

1.1 This standard prescribes the method of test for estimating parathion [*O O*-diethyl *O*-(4-nitrophenyl) phosphorothioate] foods and their milled products.

1.1.1 The method will be applicable to a minimum limit of 25 µg of parathion residues.

2. OUTLINE OF THE METHOD

2.1 The residues of parathion are extracted with benzene or isopropanol-benzene and the cleaned up solution is clarified. The nitro group of the parathion is then reduced with zinc dust and hydrochloric acid solution to amino group. The resulting water soluble compound is diazotized and coupled with *N*-(1-naphthyl) ethylene diamine to produce a magenta colour, the intensity of the colour is measured with a spectrophotometer at 555 to 560 mµ.

3. SAMPLING

3.1 The representative samples for the purpose of estimating parathion residues in various agricultural and food products shall be in accordance to the sampling procedures as prescribed in the relevant Indian Standards, wherever available.

4. EXPERIMENTAL PROCEDURE

4.1 Reagents

4.1.1 *Benzene* — analytical grade and redistilled.

4.1.2 *Adsorbent Mixture* — 10 parts of anhydrous powdered sodium sulphate, 5 parts of Atta clay, 5 parts Filter Cel, 2 parts of activated carbon.

4.1.3 *Parathion Standard Solution* — Prepare a standard solution for parathion in such a manner that 1 ml of the final solution contains 20 µg of parathion.

4.1.4 *Standard Hydrochloric Acid Solution* — 0.5 N.

4.1.5 *Sodium Nitrite Solution* — aqueous, 0.25 percent (*w/v*).

4.1.6 *Ammonium Sulphamate Solution* — aqueous, 2.5 percent (*w/v*). Dissolve 1.25 g in water and dilute to 50 ml. Prepare weekly.

4.1.7 *N*-(1-Naphthyl) *Ethylenediamine Di-Hydrochloride Solution* — aqueous, 1 percent (*w/v*). Prepare fresh solution daily and store the above in a dark coloured bottle.

4.1.8 Zinc Dust — analytical reagent-grade.

4.1.9 Isopropyl Alcohol — analytical reagent grade.

4.2 Preparation of Sample, Clean-up and Extraction of Parathion Residues

4.2.1 Firm Relatively Tough-Skinned Fruits, such as Apples, Pears, etc— Weigh 2 to 3 kg of fruit into a clean dry jar of 15 litres, so mounted that it can be turned with end-over-end tumbling action by hand-crank or motor. Add 500 ml of benzene (*see* Notes 1 and 2) and stopper with tight-fitting cork, wooden bung, or plastic screw-cap faced with gasket of sheet-cork or other suitable solvent-resisting material. Turn the jar for 5 minutes at 75 to 100 rpm (*see* Note 3). Open the jar and drain off benzene as completely as possible into 1-litre erlenmeyer flask (*see* Note 4).

NOTE 1 — In all cases; treat the benzene extract by shaking for five minutes with the adsorbent mixture, in proportion of approximately 10 g to 100 ml extract. Finally, filter through rapid, folded paper.

NOTE 2 — As benzene vapours are toxic, take precautions against inhaling them at all times.

NOTE 3 — Since apples treated for longer periods become swollen and soft, limit the time of stripping; on other firmer products time factor is not so important.

NOTE 4 — It is assumed that concentration of parathion in benzene poured off from apples is same as that in solvent retained by apples after wetting.

4.2.1.1 Aliquots of this solution may be used directly for parathion residue determination if they do not contain more than 100 mg of foreign organic matter.

4.2.1.2 In order to remove partially any colouring matter and plant waxes, add, to each 100 ml of the benzene solution, 10 g of a mixture compounded by weight by mixing 10 parts of anhydrous sodium sulphate, 5 parts Atta clay, 5 parts Filter-Cel and 2 parts Nuchar. Stopper flask, shake vigorously for five minutes, and filter on folded paper (*see* Note). Do not allow appreciable evaporation of the filtrate before determining parathion in aliquot.

NOTE — This operation dehydrates the benzene, promotes rapid filtration and colouring matter and about approximately half of the natural waxes that may be present in the products.

4.2.2 Soft Fruits, such as Plums, Tomatoes, Berries, etc— Use 1 to 2 kg sample with 300 to 500 ml of the redistilled benzene, and strip by shaking gently by hand for five minutes in a jar of suitable size.

4.2.3 Fresh, Leafy Vegetables, such as Cabbage, Lettuce, Greens, etc— Blend 100 g of the sample into 100 ml of isopropyl alcohol for two minutes in high speed blender. Add 200 ml of benzene and blend again for two minutes,

pour the mixture into a centrifuge bottle and centrifuge for 5 minutes. Transfer supernatants to a 1-litre separator. Wash blender with 50-ml portions of benzene and transfer to centrifuge bottles. Break up solids with stirring rod, stopper and shake for 2 minutes. Centrifuge as before and add solvent layer to separator. Wash and extract with equal volume of distilled water to remove *isopropyl* alcohol. Dry benzene layer with 30 g of anhydrous sodium sulphate and dilute to suitable volume with benzene and finally make the volume to 200 ml with benzene.

4.2.4 Other Products—Grind dried products, in a suitable mill and extract with benzene in a large Soxhlet-type extractor for a period of one to two hours. Satisfactory extraction may be achieved in many cases by letting ground material steep overnight in stoppered jars with measured volume of redistilled benzene. Follow the same procedure as detailed in 4.2.3.

4.3 Preparation of Standard Curve for 0-200 μ g Range—Add 2.0, 4.0, 6.0, 8.0, 10.0 ml parathion standard solution B to a series of 250-ml conical flasks. Make up the volume to 25 ml with benzene in all cases and provide a blank of 25 ml benzene. Add 20 ml of dilute hydrochloric acid and 200 mg of finely powdered zinc dust and connect to a condenser by all glass-adaptor fitted with thermometer. Place on a hot-plate at medium heat and rapidly distil off the benzene. Where the vapour temperature is about 70°C, ebullition stops, benzene is eliminated and vapour temperature falls. Do not let aqueous solution to boil at this point, disconnect the flask, add 10 ml of alcohol and reconnect to reflux condenser. Reflux for 5 minutes and cool the flask in ice water. Treat the other flasks likewise in succession.

4.3.1 Add about 100 mg of Filter Cel to each flask and filter into 50-ml volumetric flask through Whatman No. 44 or equivalent. Use long thin stirring rods during transfer of solutions and washing to filter. Take care not to exceed the volume of 45 ml in the flasks. Volume of wash water that may be used is limited to about 15 ml.

4.3.2 Add 1 ml of sodium nitrite solution to each flask, mix and let stand for 10 minutes, add 1 ml of ammonium sulphamate solution to each flask, mix, and let stand for 10 minutes. Then add 2.0 ml of reagent *N*-(1-naphthyl) ethylene diamine dihydrochloride solution, dilute to 50 ml volume and let stand for 10 minutes. Take the absorbance of the solution at 555 $m\mu$ in spectrophotometer (*see* Note). Colours should be stable for at least one hour.

NOTE—With 1-cm cell absorbance for 100 mcg of parathion is about 0.33. If available 5-cm cell may be preferred.

4.3.3 Read the absorbance in succession using blank standard as reference and plot optical density against concentration in micrograms of parathion to obtain standard curve.

4.4 Determination—Place aliquot of clarified solution containing the extract of parathion residues obtained as under **3.2** from the material under investigation in a 250-ml conical flask, the volume of aliquot taken depending upon the expected residue content. If the aliquot is less than 25 ml dilute to 25 ml with redistilled benzene, add 20 ml of 0.5 N hydrochloric acid and 200 mg of zinc dust and proceed as for the standard solution of parathion for the preparation of standard curve.

4.5 Reporting of Parathion Residues — Develop colour, read absorbance against standard blank and determine content of parathion in the aliquot taken from the standard curve and content of parathion residues in the material in terms of parts per million calculated from the same.

(Continued from page 2)

Members

SHRI V. LAKSHMINARAYANA

DR B. D. NARANG

SHRI S. N. MITRA

DR A. A. VAHORA

SHRI V. G. DESHPANDE (Alternate)

SHRI K. VISWESWARIAH

SHRI J. R. RANGASWAMY (Alternate)

DR B. L. WATTAL

SHRI V. N. BHATNAGAR (Alternate)

Representing

Directorate of Plant Protection, Quarantine & Storage
(Ministry of Food, Agriculture, Community
Development & Co-operation)

Public Analyst, Government of Haryana

Central Food Laboratory, Calcutta

CIBA of India Limited, Bombay

Central Food Technological Research Institute,
Mysore

National Institute of Communicable Diseases, Delhi

PUBLICATIONS OF INDIAN STANDARDS INSTITUTION

INDIAN STANDARDS

Over 6 500 Indian Standards covering various subjects have been issued so far. Of these, the standards belonging to the Agricultural and Food Products Group fall under the following categories:

Abattoir	Food additives
Alcoholic drinks	Food grain handling and storage
Animal feeds	Fruits and vegetables
Animal housing and equipment	Honey and by-products
Bakery and confectionery	Infant foods
Bee-keeping equipment	Laboratory animals
Beverages	Meat and meat products
Cereals, pulses and their products	Pest control equipment
Cocoa products	Pesticidal formulations
Coffee and its products	Pesticides, technical grade and general
Dairy equipment	Propagation materials
Dairy industry, layout plans	Regulated market yards
Dairy industry, methods of test	Sensory evaluation
Dairy laboratory apparatus	Spices and condiments
Dairy products	Starch derived products
Edible starch and starchy products	Sugars and by-products
Farm implements and machinery	Tea
Fish and fishery products	Tobacco products
Fish industry, sanitary conditions	Transport of live animals

OTHER PUBLICATIONS

ISI Bulletin (Published Every Month)

Single Copy	Rs 3 ⁰⁰
Annual Subscription	Rs 25 ⁰⁰

Standards: Monthly Additions

Single Copy	Rs 0 ³⁰
Annual Subscription	Rs 3 ⁰⁰

Annual Reports (from 1948-49 Onwards) ... Rs 2⁰⁰ to 5⁰⁰

Handbook of ISI Publications, 1970 ... Rs 12⁰⁰

Addendum to ' Handbook of ISI Publications, 1970 ' ... Rs 2⁵⁰

INDIAN STANDARDS INSTITUTION

Manak Bhavan, 9 Bahadur Shah Zafar Marg, New Delhi 1

Telephone : 27 01 31 (20 lines)

Telegrams : Manaksanstha

Branch Offices :

0-18 New Civil Hospital Annexe, Asarva
F Block, Unity Bldg, Narasimharaja Square
534 Sardar Vallabhbhai Patel Road
5 Chowringhee Approach
5-9-201/2-A Chirag Ali Lane
117/418 B Sarvodaya Nagar
54 General Patters Road

Ahmedabad 16
Bangalore 2
Bombay 7
Calcutta 13
Hyderabad 1
Kanpur 5
Madras 2

Telephone

—
2 76 49
35 69 44
23-08 02
5 34 35
82 72
8 72 78